

DERMATOPHYTE

REAL-TIME PCR KIT

IMPROVE YOUR DIAGNOSTICS

Rapid detection of dermatophytes in general (pan-dermatophytes) and specifically *Trichophyton rubrum* in nail samples.

SSI Diagnostica A/S
Herredsvejen 2
3400 Hillerød
Denmark
T +45 4829 9100

ssidiagnostica.com

DERMATOPHYTE

REAL-TIME PCR KIT

DESCRIPTION

The Dermatophyte Real-Time PCR Kit detects dermatophytes (pan-dermatophytes) from nail samples, specifically *T. rubrum*.

Standardized template preparation and reduced cycling allows for 2-3 hours detection time. The multiplex approach allow for simultaneous detection of several species of dermatophytes.

The kit contains all reagents needed for the analysis.

Fast identification of the problem-causing dermatophytes allows for correct treatment and thereby reduces the anti-fungal medicine prescribed plus reduces the recovery time of the patient.

PRODUCT DETAILS

Dermatophyte Real-Time PCR Kit

Article No.: 98730

Kit Size: 100 tests




Figure: Amplification curve for a *T. rubrum* positive nail sample

