


# SALMONELLA SERO-QUICK

## Serogrouping



Send the strain to a Reference Centre for serotyping.

# SEROGROUPING

with Salmonella Sero-Quick

---

## Serogroups identified

Alphabetical

O groups

---

A ~

O:2

B ~

O:4

C ~

O:7 or O:8

D ~

O:9 or O:9,46 or O:9,46,27

E ~

O:3,10 or O:1,3,19

F ~

O:11

G ~

O:13

---

If positive in D + Vi : possibly *S. Typhi*

---